



FUNCTIONAL ARMCHAIR

FREJA 239

Freja is a versatile furniture family with armchairs, sofas and chairs with different functions and sizes with tactile materials that create a varied but harmonious and uniform expression.

Freja 239 is a stepless and adjustable functional armchair with wheels and a metal frame. It offers high comfort and is easy to use. The armchair is designed for environments with high demands on cleanliness. Freja comes with removable fabric in both seat and back, as well as a gap between seat and back, which facilitates cleaning. Freja 239 is tested for 160 kg user weight according to Möbelfakta's requirements.

Freja är ritad av NC Design Studio.

ADDED VALUE

- Developed for environments with special hygiene requirements
- High seating comfort and infinitely adjustable
- Gap between seat and back for best hygiene
- Wheel
- Removable upholstery
- Optional RAL color
- Tested for public use and meets the requirements for durability & safety, set forth in the following standard, EN 16139
- Tested for 160kg

ACCESSORIES

- Metal frame in any RAL color
- Rubber fabric/hospital fabric
- Adjustable neck pillow
- Fabric-covered armrests


DIMENSIONS

Height	117 cm
Width	58 cm
Depth	79 cm
Seat height	48 cm
Seat depth	49 cm
Höjd armledare	66 cm



Any NCS-color available as an option on this product.

VIEWS



NC Nordic Care AB

Box 30 573 35 Tranås Sweden

+46(0)140-38 40 60, info@ncnordiccare.se

ncnordiccare.se



REUSE

- Removable upholstery



LEATHER

If you take care of the leather in the right way, it will become more beautiful over the years.

- Leather is heat sensitive. Therefore, direct sunlight and close contact with heating elements should be avoided.
- Vacuum the furniture to remove loose dirt and dust.
- Wipe the leather with a slightly damp cloth and a mild soapy solution can be used for cleaning.
- Water-based leather care products extend the lifespan.

METAL

- For regular cleaning, use a lightly moistened cloth with a mild detergent without solvents or abrasives.
- For heavily soiled surfaces, a cloth moistened with alcohol can be used.
- Disinfection can be done with alcohol if the surface is dried afterwards.

FABRIC

- Removable and washable fabrics are washed according to washing instructions. These are sewn with a shrink seam, wash before use for optimal fit.
- Vacuuming with a soft nozzle is often sufficient for cleaning.
- Stains should always be removed as quickly as possible. In most cases, wiping with a damp cloth is sufficient, for more difficult stains, a pH-neutral detergent or a foam detergent for professional use can be used.
- Always use white cloths and avoid wetting a woolen fabric to a greater extent, as this risks remaining stains.
- If necessary, a foam detergent can also be used for a more thorough cleaning if this is allowed in the fabric's specification.